burst kids M





Two magical hoose from hoose from

Every day is different which is why you can choose the mode that's right for you each time you brush!

Standard •


Tough on plague, gentle on gums, this mode blasts bad bacteria and other nasties with ease.

Sensitive O

Teeth and gums feeling sensitive? This mode is right for you!

Charging +

You won't see this + sign lit up very often, since a single charge can last up to four weeks! When it does light up, simply plug it in — the light will go off when it's fully charged.





Switching between modes is super easy! We'll show you how:

Press () once to turn the brush on.

Press twice to change modes.

Press a third time to turn the brush off.

